

# Cycle by Cycle Analysis

Workshop

10.08.23

By Sophia Snipes


# Schedule

|           |   |
|-----------|---|
| 2:00-3:00 | Presentation on how cycle-by-cycle analysis works |
| 3:00-3:30 | Coffee break                                      |
| 3:30-5:30 | Practical part                                    |
| 5:30-??   | Apero!  |

# Cycle-by-cycle analysis is for detecting oscillation **bursts**


Cycle-by-cycle analysis finds **bursts of oscillations** in the EEG so that they can be quantified based on their:

- Amplitudes
- Frequency
- Durations
- Spatial distribution
- Shape




# The alternatives for quantifying oscillations


Spectral power




FOOOF




Time-frequency analysis




Thresholds




# Problem 1: the EEG is non-stationary



EEG oscillation bursts  
come and go, so simple  
power is not sufficient to  
quantify changes




# Problem 2: EEG is made of both oscillations and “colored background noise”




More than white noise, colored noise makes it all too easy to find false positive oscillations in the signal.

Background also changes with brain state, independently from oscillations.



# Solution: cycle-by-cycle analysis finds bursts based on the **shape** of the signal



## Other advantages:

- High temporal resolution
- High frequency resolution
- Robust to non-sinusoidal oscillations



# What can you do with cycle-by-cycle analysis?


- Cross-frequency coupling
- Connectivity analysis
- Burst categorization by waveform
- Dissociate changes in amplitude from quantities of oscillations
- Travelling waves
- Infra-slow oscillations

**If it involves oscillations, then this is the analysis for you!**




# When NOT to use cycle-by-cycle analysis?


When interested in the background activity



When interested in isolated waves



Certain cross-frequency coupling



**How it works**

Step 1:  
**Filter** EEG  
data

Step 2:  
Detect all  
the **cycles**

Step 3:  
**Group**  
periodic  
cycles

# Step 1: Filtering EEG data

Preprocessed EEG data

Broadband filtered data

2-40 Hz

Narrowband filtered data

4-8 Hz



6-10 Hz

8-12 Hz

10-14 Hz

12-16 Hz


14-18 Hz



# Step 2: detect all the cycles

I:


Detect zero-crossings



In **narrowband** data, identify all the positive and negative zero-crossings. N.B. this is the only step on the narrowband data.

II:


Detect peaks



In the **broadband** data, detect positive and negative peaks between zero-crossings (max absolute value).

III:

Divide into cycles




Each cycle is referenced to the **negative peak**, and starts and ends with the two neighboring **positive peaks**.

# Step 3: group periodic cycles into bursts


I:

Assign properties to  
each cycle




II:

Bursts are when at least  
N cycles in a row meet  
all the criteria




# For best results, apply multiple passes of burst detection






# Criteria


Local minimum




Correct period  
(per band)




Consistent period




Consistent flank  
amplitudes



Consistent peak  
amplitudes




Monotonicity




# NOT Criteria

Minimum amplitude



The nitty-gritty  
details



# The repository

The screenshot shows the GitHub repository page for HuberSleepLab / Matcycle. The repository is public and has 0 stars, 0 forks, and 0 watchers. It has 3 branches and 2 tags. The main branch is not protected. The repository description states: "This is an adaptation of the bycycle python repository for identifying bursts in the EEG." The repository contains a README.md file, a LICENSE file, and several other files including cycle\_functions, support\_functions, .gitattributes, .gitignore, Example.m, LICENSE, PinkNoiseTest.m, README.md, and addMatcyclePaths.m. The repository was created by snipeso and has 64 commits. The repository was last updated 2 weeks ago. The repository is licensed under the MIT license. The repository has 2 releases, with the latest release being iScience, Snipes et al. 2023, published on May 20. The repository has no packages published. The repository is written in MATLAB (100.0%).

HuberSleepLab / Matcycle

Type to search

Code Issues 3 Pull requests 2 Actions Projects Wiki Security Insights Settings

Matcycle (Public)

Edit Pins Watch 0 Fork 0 Star 0

main 3 branches 2 tags Go to file Add file Code

Your main branch isn't protected

Protect this branch

snipeso Merge branch 'Lapse-causes' d293b1a 2 weeks ago 64 commits

| File               | Commit   | Time         |
|--------------------|--|--------------|
| cycle_functions    | Merge branch 'Lapse-causes'                          | 2 weeks ago  |
| support_functions  | Merge branch 'Lapse-causes'                          | 2 weeks ago  |
| .gitattributes     | Initial commit                                       | last year    |
| .gitignore         | now loops through requested repeats, and sa...       | last year    |
| Example.m          | fixes, now it works                                  | last year    |
| LICENSE            | Initial commit                                       | last year    |
| PinkNoiseTest.m    | cleanup  | 3 months ago |
| README.md          | Merge branch 'Lapse-causes'                          | 2 weeks ago  |
| addMatcyclePaths.m | copying chart, little script to add to path all s... | last year    |

README.md

## Matcycle

This is a repository of functions that detect EEG oscillation bursts based on the shape of the signal. It is a MATLAB implementation of the cycle-by-cycle analysis outlined by [Cole & Voytek, 2019](#) originally created in [python](#).

It was first used in *How and when EEG reflects changes in neuronal connectivity due to time awake*, by Snipes et al. 2023, iScience, applied in the repository [Theta\\_Bursts](#).

### How to use

About

This is an adaptation of the bycycle python repository for identifying bursts in the EEG.

Readme MIT license Activity 0 stars 0 watching 0 forks Report repository

Releases 2

iScience, Snipes et al. 2023 (Latest) on May 20

+ 1 release

Packages

No packages published Publish your first package


Languages

MATLAB 100.0%

<https://github.com/HuberSleepLab/Matcycle>

When publishing, create a "release", so you immortalize the version actually used to obtain results


# +Packages in MATLAB



- 1) Add the Matcycle folder to path
- 2) Call the main functions like so:  
`"cycy.detect_bursts()"`
- 3) Call support functions like so:  
`"cycy.utils.highpass_filter()"`



In MATLAB, if you have a folder called "private", the functions can be called by scripts within the folder above, but not by scripts outside

# The code





# Glossary



# The Inputs

## Criteria sets

```
>> CriteriaSets
```

```
CriteriaSets =
```

```
1×2 struct array with fields:
```

```
isTruePeak  
PeaksCount  
FlankConsistency  
MonotonicityInTime  
MonotonicityInAmplitude  
isProminent  
PeriodConsistency  
AmplitudeConsistency  
MinCyclesPerBurst  
VoltageNeg  
Amplitude
```

## EEG data

```
>> EEGbroadband
```

```
EEGbroadband =
```

```
struct with fields:
```

```
srate: 250  
data: [123×7501 single]  
chanlocs: [1×123 struct]
```

```
>> EEGnarrowbands
```

```
EEGnarrowbands =
```

```
1×4 struct array with fields:
```

```
srate  
data  
chanlocs
```

# The Output: Cycles

## Cycle Detection

```
>> Cycles
```

```
Cycles =
```

```
1×3381 struct array with fields:
```

```
NegPeakIdx  
PrevPosPeakIdx  
NextPosPeakIdx
```

## Cycle Properties

```
>> AugmentedCycles
```

```
AugmentedCycles =
```

```
1×3379 struct array with fields:
```

```
NegPeakIdx  
PrevPosPeakIdx  
NextPosPeakIdx  
VoltagePrevPos  
VoltageNeg  
VoltageNextPos  
isTruePeak  
PeaksCount  
PeriodPos  
PeriodNeg  
Frequency  
Amplitude  
AmplitudeRamp  
FlankConsistency  
MonotonicityInTime  
MonotonicityInAmplitude  
isProminent  
PeriodConsistency  
AmplitudeConsistency
```

# The Output: Bursts

## Bursts

```
>> Bursts
```

```
Bursts =
```

```
1×437 struct array with fields:
```

```
CyclesCount  
CycleIndexes  
NegPeakIdx  
PrevPosPeakIdx  
NextPosPeakIdx  
VoltagePrevPos  
VoltageNeg  
VoltageNextPos  
isTruePeak  
PeaksCount  
PeriodPos  
PeriodNeg  
Frequency  
Amplitude  
AmplitudeRamp  
FlankConsistency  
MonotonicityInTime  
MonotonicityInAmplitude  
isProminent  
PeriodConsistency  
AmplitudeConsistency  
Start  
End  
BurstFrequency  
Band  
ChannelIndex  
ChannelIndexLabel  
Sign  
CriteriaSetIndex
```

## Burst Cluster

```
>> BurstClusters
```

```
BurstClusters =
```


```
1×53 struct array with fields:
```

```
CyclesCount  
CycleIndexes  
NegPeakIdx  
PrevPosPeakIdx  
NextPosPeakIdx  
VoltagePrevPos  
VoltageNeg  
VoltageNextPos  
isTruePeak  
PeaksCount  
PeriodPos  
...  
ChannelIndexLabel  
Sign  
CriteriaSetIndex  
ClusterBurstsIdx  
ClusterChannelIndexes  
ClusterChannelLabels  
ClusterStarts  
ClusterEnds  
ClusterCycleCounts  
ClusterSigns  
ClusterFrequency  
ClusterAmplitude  
ClusterAmplitudeSum  
ClusterPeaks  
ClusterStart  
ClusterEnd  
ClusterGlobality
```

Cycle information  
of reference burst



# The process






# How to choose criteria and thresholds?


Does it correctly capture obvious bursts?

Why did it miss this burst?



Does it correctly reject background activity?



Was a criteria I didn't use better?



# Diagnostics: which criteria were more aggressive?




If one criteria is solely responsible for disqualifying a lot of cycles, threshold might be too strict

Could have omitted



# Diagnostics: properties distributions



Use to decide  
where to set  
thresholds

# Practical Part

# Recommended Steps

1. Download the repository:  
<https://github.com/HuberSleepLab/Matcycle>
2. Run **SingleBandExample.m**. To understand how the analysis works, try adapting:
  - The EEG channel (if you have it, try with your own data!)
  - The narrow band frequency range
  - Criteria set thresholds
3. Run & adapt **MultiChannelExample.m** for your own data
4. Exercise for home: how long does it take to fall asleep?
  - Automatically quantify the interval between the last alpha burst and the first sleep spindle for your data

# Bibliography

- Donoghue, T., Haller, M., Peterson, E. J., Varma, P., Sebastian, P., Gao, R., ... & Voytek, B. (2020). Parameterizing neural power spectra into periodic and aperiodic components. *Nature neuroscience*, 23(12), 1655-1665.
- Cole, S., & Voytek, B. (2019). Cycle-by-cycle analysis of neural oscillations. *Journal of neurophysiology*, 122(2), 849-861.
- Snipes, S., Meier, E., Meissner, S. N., Landolt, H. P., & Huber, R. (2023). How and when EEG reflects changes in neuronal connectivity due to time awake. *iScience*.